Specifications

Туре	Class I , Type B equipment	
Photographic field	0.25 × 0.55mm	
Photographic capability	Automatic	
Photographic location	Center,Peripheral locations(0,45,135,180,225,315 degrees.)	
Central corneal thickness	Measurable	
Illumination	LED	
Analyzing method	Auto Center Method / Auto F Center Method / Auto Trace Method / Center Method / F Center Method / Trace Method	
Analysis data	Cell density (mm ²), Coefficient of variation, Standard deviation, Percent of Hexagonality, Number of analyzed cells, Average cell area (μ m ²), Maximum cell area (μ m ²), Minimum cell area (μ m ²), Distribution of number of sides histogram (%), Distribution of area histogram (%)	
Integrated monitor	10.1" Capacitive touch panel with 180 degrees flexible horizontal & vertical movement capability	
External Interface	EMR connection / DICOM capability	
External ports	USB3.0 x 4 Type A , RJ-45	
Power	AC100-240V , 50 / 60Hz	
Power consumption	100VA	
Dimensions	$310 (W) \times 459 (D) \times 451 (H) mm ($ with the monitor at the rear-side $)$	
Weight	19.6kg	



%Specifications subject to change without notice.

Product Name : Konan Specular Microscope XX



CellChek 20-1

0000.23.11B



Specular Microscopy & Optical Pachymetry



0000

[Product Name : Konan Specular Microscope XX]

Equipped with "Auto Center Method" and "Auto F Center Method" as well as "Auto Trace Method" as standard!

Image Capturing

Enhanced Image Capturing Capability

By comparison with our conventional model*, the capturing time shortened by approx. 11%, the image field widened by approx. 43% to 0.25mm x 0.55mm and the resolution improved by approx. 25%. Approx. 35% shorter booting time enables even faster exam. *Compared with CellChekSL

Automated Capturing Retry Feature

CellChek 20-1 automatically retries the preset number of the captures if it is unable to capture cell images due to blinking or eyelashes.

Simultaneous Non-Contact Pachymetry Central corneal thickness can be measured simultaneously.

Captured Location Confirmation

The location captured can be confirmed on the anterior segment image. It can be used for pre and post surgery assessment.



One central and 6 peripheral locations can be captured.

In normal capturing mode, multiple corneal endothelial cell images can be captured and analyzed in a single examination.



Sample screen shot of the normal capturing mode

Manual Capturing

It is useful for abnormal corneas that require extra effort for imaging.



Simple to use fully automated OD/OS endothelium capturing, analysis, printing and exporting, with a touch of screen.



or

 $\langle \rangle$

Export

Complete

Capture result	
	OD
\gg	
AUTO	2

Touch capturing button Alignment or center of the pupil. & Focus

AUTO







OD/OS results



 \ll

AUTO

OD/OS analysis





AUTO

endothelium capturing





OS



Auto Center Method, Auto F Center Method, Auto Trace Method, Center Method, F Center Method, Trace Method

The "Center Method" which has been widely used in clinical practice and clinical trials around the world for many years, and the "F Center Method", are automated. "Auto Center Method" and "Auto F Center Method" are available. The required analysis results can be obtained with a touch of the screen. "Auto Trace Method" is also available, and can be used in different ways depending on the situation.





Sample image of Auto Center Method

Sample image of Auto F Center Method



Analysis result screen

Output & Data Management

Output

Compatible with various types of external printers. DICOM capability.

Filing Features

Captured image and analysis data are automatically saved. Data retrieval can be performed as needed. Data transfer from our conventional model is possible.

Patient Progression

The patient progression of patient data can be displayed.

Installation Layout

Flexible Monitor Direction

The large 10.1-inch touch monitor can be freely moved to various orientations. It can be moved in accordance with the operating environment to provide support to the examinee from various angles.

The main unit is compact, lightweight, and has no moving parts during measurement, allowing installation layout according to the space requirements.













Sample image of the patient progression



*The expressions in this brochure are based on comparison with our previous product, CellChek SL.